

DMR plus Specifications

CR Ready

Generator

Constant potential high frequency
kV range: 22 - 49kV
mAs range: 3 - 600mAs
Input line: 200 - 415 VAC, 30 amps

X-ray tube

Dual track, Rotating anode
Target type: Molybdenum / Rhodium
Speed: 9000 RPM max
Focal spots: 0.1 x 0.3 mm
Filtration: 0.03mm Molybdenum / 0.025mm Rhodium / 1.0mm Aluminum

C-arm

Automatic dual compression
End of exposure automatic decompression
Manual compression
Compression display readout, force and thickness
Angulation display readout
Motorized vertical positioning
Rotation + / - 180°
Automatic collimation
SID: variable 560 to 660 mm
Stereotactic: upgradeable

Operator control console

Operator modes: AOP (Automatic optimization of parameters)
Auto-time and manual
Operator shield

Standard accessories

Dual buckys 18 x 24 cm and 24 x 30 cm
Grids: linear
Dual compression paddles 18 x 24 cm and 24 x 30 cm
Spot compression paddle
Magnification platform (1.5 and 1.8)
ID flasher
Face shield
Accessory cart
Manuals: operator, service and schematics

Optional accessories

Localization paddle and crosshair
Stereotactic: positioner and controller

