

# Performa Specifications

CR Ready

## Generator

Constant potential high frequency  
kV range: 23- 35kV  
mAs range: 4 - 500mAs  
Input line: 220 - 240 VAC, 30 amps, 50/60Hz

## X-ray tube

Rotating anode  
Target type: Molybdenum  
Focal spots: 0.1 x 0.3 mm  
Filtration: 0.03 mm Molybdenum / 0.025 mm Rhodium

## C-arm

Automatic dual compression  
End of exposure automatic decompression  
Manual compression  
Compression display readout, force and thickness  
Angulation display readout  
Motorized vertical positioning  
Manual rotation: +135 and -180°

## Operator control console

Operator modes: full auto, semi auto and manual  
Operator shield  
ID Flasher

## Standard accessories

Dual bucky: 18 x 24 cm and 24 x 30 cm  
Grids: linear  
Dual compression paddles 18 x 24 cm and 24 x 30 cm  
Spot contact compression paddle  
Magnification platform (1.6 and 1.8)  
Spot magnification compression paddle  
Face shield  
Manuals: operator, service and schematics

## Optional accessories

Accessory cabinet  
Localization paddle  
ROC Grids

